

MiniScan

Radio-TLC and Radio-HPLC-System

MiniScan

MiniScan is a versatile, low-cost radio-TLC scanner for the reliable detection of radioisotopes on narrow strips and plates.



Technology

MiniScan is a versatile TLC scanner for the reliable detection of radioisotopes on narrow strips and plates. The system is ideal for routine quality control of [^{18}F]FDG, $^{99\text{m}}\text{Tc}$ and ^{123}I radiopharmaceuticals.

A complete MiniScan system consists of a moving stage, a FlowCount system and a PMT based detector. FlowCount is a radioisotope HPLC detection system, which is compatible with all HPLC systems. MiniScan uses various interchangeable NaI photomultiplier detectors for measurements of most isotopes including ^{18}F , ^{125}I , ^{131}I , $^{99\text{m}}\text{Tc}$, and ^{111}In . Several scan speeds and variable detector slits allow the detectors to measure a wide range of activities from 10 nCi to 100 μCi . Analog and digital signals are provided for interfacing MiniScan with existing chromatography data systems. *RaPET Chromatography Software* can also be used for your data collection and report generation requirements.

Applications

- Routine quality control of [^{18}F]FDG, $^{99\text{m}}\text{Tc}$, ^{123}I radiopharmaceuticals
- TLC of radiopharmaceuticals labeled with gamma, beta, and alpha emitters
- ^{14}C verification in pharmaceutical and toxicology studies
- In-process TLC analysis of reaction mixtures

RaPET Chromatography Software

RaPET Chromatography Software is a comprehensive chromatography data collection and analysis package. This evaluation software is very reliable and easy-to-use. It furthermore consists of a GMP database ensuring GMP compliant documentation and adhering to 21 CFR part 11.

Features and benefits

- Easy-to-use system with maximum flexibility
- Easily adaptable for dual use as an HPLC detection system: FlowCount (requires flow cells and detector holder)
- Compatible with all chromatography data systems and *RaPET Chromatography Software*
- Range of detectors to suit many applications
- Capable of monitoring two detectors simultaneously (Radio-HPLC mode)
- Easy system setup and maintenance

Models

MS-1000F

MiniScan TLC Radiochromatography system with single PMT FlowCount base unit and scanning stage. Variable scan speeds, collimated detector holder, analog and digital outputs.

MS-2000F

MiniScan TLC and HPLC Radiochromatography system with FlowCount dual detector unit and scanning stage. Variable scan speeds, collimated detector holder, analog and digital outputs.

MS-2000FP

MiniScan TLC and HPLC Radiochromatography system with FlowCount dual PMT base unit and scanning stage. Variable scan speeds, collimated detector holder, analog and digital outputs.

MiniScan


Detectors

Two types of detectors are available for the MiniScan system, a standard Na/I crystal and a plastic scintillator. These interchangeable detectors provide the system with the flexibility to measure a wide range of isotopes and activities. See the FlowCount documentation for more information about radio-HPLC detectors.

- The FC-3100 Na/I PMT based detector is a low energy gamma (10 - 60 keV) detector used primarily for ^{125}I .
- The FC-3200 Na/I PMT based detector is a high energy gamma (>60 keV) detector used in most nuclear medicine applications.
- The FC-3600 Plastic Scintillator/PMT based detector is ideal for the detection of ^{32}P , ^{90}Y and other high energy beta emitters.

Upgrade your MiniScan to a Radio-HPLC detection system

Check our FlowCount flyer for Radio-HPLC detectors, detector holders and flow cells.



Eckert & Ziegler Radiopharma, Inc.

25 Upton Drive
Wilmington, MA 01887
USA

Phone: + 1 508 497 0060

Fax: + 1 508 497 0061

info@radiopharma-inc.com

www.radiopharma.com